

Product Name: Rel B Rabbit Monoclonal Antibody
Catalog #: AMRe04107



Summary

Production Name	Rel B Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purification

Immunogen

Gene Name	RELB
Alternative Names	RELB; Transcription factor RelB; I-Rel
Gene ID	5971
SwissProt ID	Q01201.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 62 kDa; Observed MW: 62 kDa

Background

Product Name: Rel B Rabbit Monoclonal Antibody
Catalog #: AMRe04107




Transcription factors of the nuclear factor κ B (NF- κ B)/Rel family play a pivotal role in inflammatory and immune responses. Involved in many biological processes such as inflammation, immunity, differentiation, cell growth, tumorigenesis and apoptosis.


Research Area

Cell Biology

Image Data



Western blot analysis of Rel B in 3T3, C6 lysates using Rel B antibody.



Western blot analysis of Rel B in Raji lysates using Rel B antibody.

Note

For research use only.