

**Product Name: Phospho-mTOR (Ser2481) (10P6) Rabbit Monoclonal Antibody**  
**Catalog #: AMRe05952**

---

## Summary

|                        |                                                          |
|------------------------|----------------------------------------------------------|
| <b>Production Name</b> | Phospho-mTOR (Ser2481) (10P6) Rabbit Monoclonal Antibody |
| <b>Description</b>     | Rabbit Monoclonal Antibody                               |
| <b>Host</b>            | Rabbit                                                   |
| <b>Application</b>     | WB                                                       |
| <b>Reactivity</b>      | Human,Mouse,Rat                                          |

## Performance

|                     |                                                                                                               |
|---------------------|---------------------------------------------------------------------------------------------------------------|
| <b>Conjugation</b>  | Unconjugated                                                                                                  |
| <b>Modification</b> | Phospho Antibody                                                                                              |
| <b>Isotype</b>      | IgG                                                                                                           |
| <b>Clonality</b>    | Monoclonal                                                                                                    |
| <b>Form</b>         | Liquid                                                                                                        |
| <b>Storage</b>      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.                      |
| <b>Buffer</b>       | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. |
| <b>Purification</b> | Affinity purification                                                                                         |

## Immunogen

|                          |         |
|--------------------------|---------|
| <b>Gene Name</b>         | MTOR    |
| <b>Alternative Names</b> |         |
| <b>Gene ID</b>           | 2475.0  |
| <b>SwissProt ID</b>      | P42345. |

## Application

|                         |                  |
|-------------------------|------------------|
| <b>Dilution Ratio</b>   | WB 1:1000-1:5000 |
| <b>Molecular Weight</b> | 289 kDa          |

## Background


---

**Product Name: Phospho-mTOR (Ser2481) (10P6) Rabbit  
Monoclonal Antibody  
Catalog #: AMRe05952**

---

## Research Area

## Image Data



Western blot analysis of Phospho-mTOR (Ser2481) in C6 lysates using Phospho-mTOR (Ser2481) antibody.

## Note

For research use only.