

Product Name: BCL2L12 (18T13) Rabbit Monoclonal Antibody
Catalog #: AMRe07503



Summary

Production Name	BCL2L12 (18T13) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,IHC-P,ICC/IF,FC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	BCL2L12
Alternative Names	BCL2 like 12; BCL2L12; BPR;
Gene ID	83596.0
SwissProt ID	Q9HB09.

Application

Dilution Ratio	WB 1:2000-1:20000, IHC-P/IF-P 1:200-1:1000, ICC/IF , FCM 1:20-1:100
Molecular Weight	37kDa

Product Name: BCL2L12 (18T13) Rabbit Monoclonal Antibody
Catalog #: AMRe07503




Background

Belongs to the Bcl-2 family.

Research Area

Image Data



Western blot analysis of BCL2L12 expression in MCF7 cell lysate.

Note

For research use only.