

Product Name: Desmoglein 3 (7G9) Rabbit Monoclonal Antibody
Catalog #: AMRe09926



Summary

Production Name	Desmoglein 3 (7G9) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	IHC-P,IF-P
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	DSG3
Alternative Names	CDHF6; Desmoglein3; DSG3; PVA;
Gene ID	1830.0
SwissProt ID	P32926.

Application

Dilution Ratio	IHC-P/IF-P 1:100-1:200
Molecular Weight	108kDa

Product Name: Desmoglein 3 (7G9) Rabbit Monoclonal Antibody
Catalog #: AMRe09926




Background

Component of intercellular desmosome junctions. Involved in the interaction of plaque proteins and intermediate filaments mediating cell-cell adhesion. Component of intercellular desmosome junctions. Involved in the interaction of plaque proteins and intermediate filaments mediating cell-cell adhesion.

Research Area

Image Data



Western blot analysis of Desmoglein 3 expression in A431 cell lysate.

Note

For research use only.