

Product Name: FUBP1 (9R5) Rabbit Monoclonal Antibody
Catalog #: AMRe11179

Summary

Production Name	FUBP1 (9R5) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,IHC-P,ICC/IF,FC,IP,IF-P
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	FUBP1
Alternative Names	
Gene ID	8880.0
SwissProt ID	Q96AE4. A synthetic peptide of human FUBP1

Application


Dilution Ratio	WB 1:1000, IHC-P/IF-P 1:1000-1:10000, ICC/IF 1:200-1:500, FCM 1:200-1:500, IP 1:20-1:50
Molecular Weight	68kDa

Product Name: FUBP1 (9R5) Rabbit Monoclonal Antibody
Catalog #: AMRe11179

Background

Research Area

Image Data



Western blot detection of FUBP1/GBP in K562,Rat Brain,HeLa cell lysates using FUBP1/GBP antibody(1:1000 diluted).

Note

For research use only.