

Product Name: GM130 (6E7) Rabbit Monoclonal Antibody
Catalog #: AMRe11527

Summary

Production Name	GM130 (6E7) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,IHC-P,ICC/IF,FC,IP,IF-P
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol.
Purification	Affinity purification

Immunogen

Gene Name	GOLGA2
Alternative Names	
Gene ID	2801.0
SwissProt ID	Q08379.

Application


Dilution Ratio	WB 1:1000-1:5000, IHC-P/IF-P 1:100-1:500, ICC/IF 1:50-1:200, FCM 1:20-1:50, IP 1:20-1:50
Molecular Weight	113 kDa

Product Name: GM130 (6E7) Rabbit Monoclonal Antibody
Catalog #: AMRe11527

Background

Research Area

Image Data



Western blot analysis of GM130 in HeLa, lysates using GM130 antibody.

Note

For research use only.