

Product Name: Histone H1.0 (10F13) Rabbit Monoclonal Antibody
Catalog #: AMRe12052

Summary

Production Name	Histone H1.0 (10F13) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,IHC-P,ICC/IF,IF-P
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	H1-0
Alternative Names	
Gene ID	3005.0
SwissProt ID	P07305. A synthetic peptide of human Histone H1.0

Application


Dilution Ratio	WB 1:1000-1:5000, IHC-P/IF-P 1:100-1:200, ICC/IF 1:100-1:200
Molecular Weight	21kDa

Background

Product Name: Histone H1.0 (10F13) Rabbit Monoclonal Antibody
Catalog #: AMRe12052

Research Area

Image Data



Western blot detection of Histone H1.0 in C6,CHO-K1 cell lysates using Histone H1.0 antibody(1:1000 diluted).

Note

For research use only.