

**Product Name: IKK alpha/beta (10K2) Rabbit
Monoclonal Antibody
Catalog #: AMRe12472**

Summary

Production Name	IKK alpha/beta (10K2) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol.
Purification	Affinity purification

Immunogen

Gene Name	CHUK/IKBKB
Alternative Names	
Gene ID	1147/3551
SwissProt ID	O15111/O14920.

Application


Dilution Ratio	WB 1:1000, IP 1:20-1:50
Molecular Weight	85 kDa

Background

**Product Name: IKK alpha/beta (10K2) Rabbit
Monoclonal Antibody
Catalog #: AMRe12472**

Research Area

Image Data



Western blot analysis of IKK alpha/beta in 293T lysates using IKK alpha/beta antibody.

Note

For research use only.