

Product Name: LDHA (4H19) Rabbit Monoclonal Antibody
Catalog #: AMRe13262

Summary

Production Name	LDHA (4H19) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,IHC-P,ICC/IF,FC,IF-P
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type
Buffer	preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	LDHA
Alternative Names	LDHA; LDH-M; LDH1; PIG19;
Gene ID	3939.0
SwissProt ID	P00338.

Application

Dilution Ratio	WB 1:1000-1:5000, IHC-P/IF-P 1:50-1:200, ICC/IF 1:20-1:50, FCM 1:20
Molecular Weight	37kDa

Product Name: LDHA (4H19) Rabbit Monoclonal Antibody
Catalog #: AMRe13262




Background

Fermentation; pyruvate fermentation to lactate; (S)-lactate from pyruvate: step 1/1.

Research Area

Image Data



Western blot analysis of LDHA expression in MCF7 cell lysate.

Note

For research use only.