

**Product Name: MSI2 (13U1) Rabbit Monoclonal Antibody**  
**Catalog #: AMRe14176**

---

## Summary

<b>Production Name</b>	MSI2 (13U1) Rabbit Monoclonal Antibody
<b>Description</b>	Rabbit Monoclonal Antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IHC-P,ICC/IF,FC,IP
<b>Reactivity</b>	Human,Mouse,Rat

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type
<b>Buffer</b>	preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
<b>Purification</b>	Affinity purification

## Immunogen

<b>Gene Name</b>	MSI2
<b>Alternative Names</b>	MSI2H; MGC3245; MSI2;
<b>Gene ID</b>	124540.0
<b>SwissProt ID</b>	Q96DH6.

## Application

<b>Dilution Ratio</b>	WB 1:1000, IHC-P/IF-P 1:200-1:500, ICC/IF 1:100, FCM 1:100, IP 1:50
<b>Molecular Weight</b>	35kDa

**Product Name: MSI2 (13U1) Rabbit Monoclonal Antibody**  
**Catalog #: AMRe14176**


---

## Background

Msi2 (musashi homolog 2), also known as MSI2H, is a 328 amino acid protein that localizes to the cytoplasm and contains two RRM (RNA recognition motif) domains. Expressed ubiquitously at low levels, Msi2 functions as an RNA binding protein that, by regulating the expression of target mRNAs, is thought to play a role in the proliferation and maintenance of stem cells within the central nervous system. Msi2 is subject to posttranslational phosphorylation and is upregulated in response to brain injury, suggesting a role in healing and brain tissue regeneration. RNA binding protein that regulates the expression of target mRNAs at the translation level. May play a role in the proliferation and maintenance of stem cells in the central nervous system (By similarity).

## Research Area

## Image Data



Western blot analysis of MSI2 expression in T47 D cell lysate.

## Note

For research use only.