

Product Name: Occludin (10Z10) Rabbit Monoclonal Antibody
Catalog #: AMRe15092



Summary

Production Name	Occludin (10Z10) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,IHC-P,ICC/IF,IP,FC,IF-P
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type
Buffer	preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	OCLN
Alternative Names	OCLN; Occludin; BLCPMG;
Gene ID	100506658.0
SwissProt ID	Q16625.

Application

Dilution Ratio	WB 1:1000, IHC-P/IF-P 1:200, ICC/IF 1:100, IP 1:20, FCM 1:50
Molecular Weight	59kDa

Product Name: Occludin (10Z10) Rabbit Monoclonal Antibody
Catalog #: AMRe15092




Background

May play a role in the formation and regulation of the tight junction (TJ) paracellular permeability barrier. It is able to induce adhesion when expressed in cells lacking tight junctions. May play a role in the formation and regulation of the tight junction (TJ) paracellular permeability barrier. It is able to induce adhesion when expressed in cells lacking tight junctions.

Research Area

Image Data



Western blot analysis of Occludin expression in HepG2 cell lysate.

Note

For research use only.