

**Product Name: Serum Amyloid A (1K10) Rabbit
Monoclonal Antibody
Catalog #: AMRe17770**

Summary

Production Name	Serum Amyloid A (1K10) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,IHC-P,IF-P
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type
Buffer	preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	SAA1
Alternative Names	Amyloid fibril protein AA; Amyloid fibrils; Amyloid protein A; Fibrils; PIG4; SAA 1; SAA 2; Serum amyloid A protein precursor;
Gene ID	6288.0
SwissProt ID	P0DJ18.

Application

Dilution Ratio	WB 1:1000-1:5000, IHC-P/IF-P 1:100-1:200
Molecular Weight	14kDa

**Product Name: Serum Amyloid A (1K10) Rabbit
Monoclonal Antibody**
Catalog #: AMRe17770




Background

Major acute phase reactant. Apolipoprotein of the HDL complex. Major acute phase protein.

Research Area

Image Data



Western blot analysis of Serum Amyloid A expression in human plasma lysate.

Note

For research use only.