

Product Name: TNFSF9 (13Z2) Rabbit Monoclonal Antibody
Catalog #: AMRe19094

Summary

Production Name	TNFSF9 (13Z2) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type
Buffer	preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	TNFSF9
Alternative Names	4 1BB L; 4 1BBL; 4-1BB ligand; 4-1BBL; Cd137l; Cd157l; Ly63l; Tnfsf9; TNLG5A;
Gene ID	8744.0
SwissProt ID	P41273.

Application

Dilution Ratio	WB 1:1000, IP 1:50
Molecular Weight	27kDa

Product Name: TNFSF9 (13Z2) Rabbit Monoclonal Antibody
Catalog #: AMRe19094




Background

Cytokine that binds to TNFRSF9. Induces the proliferation of activated peripheral blood T-cells. May have a role in activation-induced cell death (AICD). May play a role in cognate interactions between T-cells and B-cells/macrophages. Cytokine that binds to TNFRSF9. Induces the proliferation of activated peripheral blood T-cells. May have a role in activation- induced cell death (AICD). May play a role in cognate interactions between T-cells and B-cells/macrophages.

Research Area

Image Data



Western blot analysis of TNFSF9 expression in HEK293 cell lysate.

Note

For research use only.