

Product Name: TFEB Rabbit Monoclonal Antibody
Catalog #: AMRe21127



Summary

Production Name	TFEB Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,IHC,IF,IP,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG,Kappa
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Purification	Protein A

Immunogen

Gene Name	TFEB
Alternative Names	BHLHE35
Gene ID	7942.0
SwissProt ID	P19484.

Application

Dilution Ratio	IHC 1:200-1:1000;WB 1:2000-1:10000;IF 1:200-1:1000;ELISA 1:5000-1:20000;IP 1:50-1:200;
Molecular Weight	Calculated MW:53kD;Observed MW:65-70kD

Background


Cell localization:Cytoplasm, Nucleus.

Product Name: TFEB Rabbit Monoclonal Antibody
Catalog #: AMRe21127



Research Area

Image Data



Rat thymus whole cell lysates were separated by 10% SDS-PAGE, and the membrane was blotted with primary antibody(1:1000). The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody.

Note

For research use only.