

**Product Name: Acetyl Coenzyme A Carboxylase beta  
Rabbit Monoclonal Antibody  
Catalog #: AMRe21271**

---



## Summary

<b>Production Name</b>	Acetyl Coenzyme A Carboxylase beta Rabbit Monoclonal Antibody
<b>Description</b>	Rabbit Monoclonal Antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IHC,IF,IP,ELISA
<b>Reactivity</b>	Human,Mouse,Rat

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG,Kappa
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
<b>Purification</b>	Protein A

## Immunogen

<b>Gene Name</b>	ACACB
<b>Alternative Names</b>	ACACB,Acetyl CoA carboxylase 2,ACC beta,ACC2,ACCB,AcetylCoA carboxylase 2,ACCBeta,ACCβ,ACC β,
<b>Gene ID</b>	32.0
<b>SwissProt ID</b>	O00763.

## Application

<b>Dilution Ratio</b>	IHC 1:1000-1:4000;WB 1:2000-1:10000;IF 1:200-1:1000;ELISA 1:5000-1:20000;IP 1:50-1:200;
<b>Molecular Weight</b>	Calculated MW:227kD;Observed MW:265kD

**Product Name: Acetyl Coenzyme A Carboxylase beta  
Rabbit Monoclonal Antibody  
Catalog #: AMRe21271**

---




## Background

Cell localization: Mitochondrion. Acetyl-CoA carboxylase (ACC) is a complex multifunctional enzyme system. ACC is a biotin-containing enzyme which catalyzes the carboxylation of acetyl-CoA to malonyl-CoA, the rate-limiting step in fatty acid synthesis. ACC-beta is thought to control fatty acid oxidation by means of the ability of malonyl-CoA to inhibit carnitine-palmitoyl-CoA transferase I, the rate-limiting step in fatty acid uptake and oxidation by mitochondria. ACC-beta may be involved in the regulation of fatty acid oxidation, rather than fatty acid biosynthesis. There is evidence for the presence of two ACC-beta isoforms. [provided by RefSeq, Jul 2008],

## Research Area

## Image Data



PC-12 whole cell lysates were separated by 10% SDS-PAGE, and the membrane was blotted with primary antibody(1:1000).  
The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody.

## Note

For research use only.