

Product Name: MKLP1 Rabbit Polyclonal Antibody
Catalog #: APRab00001



Summary

Production Name	MKLP1 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,ICC/IF
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification	Affinity Chromatography

Immunogen

Gene Name	KIF23
Alternative Names	CHO1; KNSL5; MKLP1; MKLP-1
Gene ID	9493
SwissProt ID	Q02241.

Application

Dilution Ratio	WB: 1:500-1:1000 IF: 1:50-1:200
Molecular Weight	Calculated MW: 110 kDa; Observed MW: 98 kDa

Background

Plus-end-directed motor enzyme that moves antiparallel microtubules in vitro. Localizes to the interzone of mitotic

Product Name: MKLP1 Rabbit Polyclonal Antibody
Catalog #: APRab00001




spindles.

Research Area

Cell Biology

Image Data



Western blot analysis of MKLP1 in HeLa lysates using MKLP1 antibody.

Note

For research use only.