

Product Name: PDGFR beta Rabbit Polyclonal Antibody
Catalog #: APRab00010



Summary

Production Name	PDGFR beta Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,IHC-P,ICC/IF,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification	Affinity Chromatography

Immunogen

Gene Name	PDGFRB
Alternative Names	PDGFRB; PDGFR; PDGFR1; Platelet-derived growth factor receptor beta; PDGF-R-beta; PDGFR-beta; Beta platelet-derived growth factor receptor; Beta-type platelet-derived growth factor receptor; CD140 antigen-like family member B; Platelet-deri
Gene ID	5159
SwissProt ID	P09619.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IF: 1:50-1:200 IP: 1:20
Molecular Weight	Calculated MW: 124 kDa; Observed MW: 170 kDa

Product Name: PDGFR beta Rabbit Polyclonal Antibody
Catalog #: APRab00010




Background

This gene encodes a cell surface tyrosine kinase receptor for members of the platelet-derived growth factor family. These growth factors are mitogens for cells of mesenchymal origin. The identity of the growth factor bound to a receptor monomer determines whether the functional receptor is a homodimer or a heterodimer, composed of both platelet-derived growth factor receptor alpha and beta polypeptides.


Research Area

Cardiovascular

Image Data




Western blot analysis of PDGF Receptor beta in SH-SY5Y lysates using PDGFR beta antibody.



Immunohistochemistry analysis of paraffin-embedded Human uterus using PDGF Receptor beta antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Product Name: PDGFR beta Rabbit Polyclonal Antibody
Catalog #: APRab00010



Immunofluorescence analysis of PDGFR beta in 3T3 using PDGF Receptor beta antibody.

Note

For research use only.