

Product Name: Calretinin Rabbit Polyclonal Antibody
Catalog #: APRab00206



Summary

Production Name	Calretinin Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,IHC-P,ICC/IF,FC
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification	Affinity Chromatography

Immunogen

Gene Name	CALB2
Alternative Names	Calretinin; CR; 29 kDa calbindin
Gene ID	794
SwissProt ID	P22676.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IF: 1:50-1:200 FC: 1:50-1:100
Molecular Weight	Calculated MW: 32 kDa; Observed MW: 29 kDa

Background

Calretinin is a calcium-binding protein which is abundant in auditory neurons.


Product Name: Calretinin Rabbit Polyclonal Antibody
Catalog #: APRab00206



Research Area

Tags & Cell Markers

Image Data



Western blot analysis of Calretinin in SH-SY5Y lysates using Calretinin antibody.

Note

For research use only.