

Product Name: TCTP Rabbit Polyclonal Antibody
Catalog #: APRab00222



Summary

Production Name	TCTP Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,IHC-P,ICC/IF,FC,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification	Affinity Chromatography

Immunogen

Gene Name	TPT1
Alternative Names	Translationally-controlled tumor protein; TCTP; Fortilin; Histamine-releasing factor; HRF; p23
Gene ID	7178
SwissProt ID	P13693.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IF: 1:50-1:200 IP: 1:20 FC: 1:50-1:100
Molecular Weight	Calculated MW: 20 kDa; Observed MW: 20 kDa

Background

Product Name: TCTP Rabbit Polyclonal Antibody
Catalog #: APRab00222




Involved in calcium binding and microtubule stabilization.

Research Area

Immunology

Image Data



Western blot analysis of TCTP Receptor in HeLa lysates using TCTP antibody.

Note

For research use only.