

Product Name: GST3 Rabbit Polyclonal Antibody
Catalog #: APRab00343



Summary

Production Name	GST3 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,IHC-P,ICC/IF
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification	Affinity Chromatography

Immunogen

Gene Name	GSTP1
Alternative Names	GSTP1; FAEES3; GST3; Glutathione S-transferase P; GST class-pi; GSTP1-1
Gene ID	2950
SwissProt ID	P09211.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IF: 1:50-1:200
Molecular Weight	Calculated MW: 23 kDa; Observed MW: 23 kDa

Background

Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles.


Product Name: GST3 Rabbit Polyclonal Antibody
Catalog #: APRab00343



Research Area

Cardiovascular

Image Data



Western blot analysis of GST3 / GST pi in PC-3 lysates using GST3 antibody.

Note

For research use only.