

Product Name: STAT6 Rabbit Polyclonal Antibody
Catalog #: APRab00349



Summary

Production Name	STAT6 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,IHC-F,IHC-P,ICC/IF,FC,IP
Reactivity	Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purification

Immunogen

Gene Name	STAT6
Alternative Names	STAT6; Signal transducer and activator of transcription 6; IL-4 Stat
Gene ID	6778
SwissProt ID	P42226.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IF: 1:50-1:200 IP: 1:20 FC: 1:50-1:100
Molecular Weight	Calculated MW: 94 kDa; Observed MW: 110 kDa

Background

STAT6 transcription factor of the STAT family. Plays a central role in IL4-mediated biological responses. Induces the

Product Name: STAT6 Rabbit Polyclonal Antibody
Catalog #: APRab00349




expression of BCL2L1/BCL-X(L), which is responsible for the anti-apoptotic activity of IL4.

Research Area

Signal Transduction

Image Data



Western blot analysis of STAT6 in Raw264.7 lysates using STAT6 antibody.

Note

For research use only.