

Product Name: kappa Actin Rabbit Polyclonal Antibody
Catalog #: APRab00393



Summary

Production Name	kappa Actin Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Chromatography

Immunogen

Gene Name	ACTBL2
Alternative Names	ACT
Gene ID	345651
SwissProt ID	Q562R1.

Application

Dilution Ratio	WB: 1:500-1:1000 ELISA: 1:10000
Molecular Weight	Calculated MW: 42 kDa; Observed MW: 42 kDa

Background

Actins are highly conserved proteins that are involved in various types of cell motility and are ubiquitously expressed in all

Product Name: kappa Actin Rabbit Polyclonal Antibody
Catalog #: APRab00393




eukaryotic cells.


Research Area

Tags & Cell Markers

Image Data



Western blot analysis of kappa Actin in various lysates using kappa Actin antibody.



Western blot analysis of kappa Actin in the HUVEC lysates using ACTBL2 antibody.

Note

For research use only.