

Product Name: PEPC Rabbit Polyclonal Antibody
Catalog #: APRab00439



Summary

Production Name	PEPC Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Plants

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	PGC
Alternative Names	PGC; Gastricsin; Pepsinogen C
Gene ID	5225
SwissProt ID	P20142.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 42 kDa; Observed MW: 105-109 kDa

Background

Hydrolyzes a variety of proteins.


Product Name: PEPC Rabbit Polyclonal Antibody
Catalog #: APRab00439



Research Area

Cardiovascular

Image Data



Western blot analysis of PEPC in Maize lysates using PEPC antibody.

Note

For research use only.