

Product Name: Toll-Like Receptor 4 Rabbit Polyclonal Antibody
Catalog #: APRab00443



Summary

Production Name	Toll-Like Receptor 4 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	IHC-P
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	TLR4
Alternative Names	TLR4; Toll-like receptor 4; hToll; CD antigen CD284
Gene ID	7099
SwissProt ID	O00206.

Application

Dilution Ratio	IHC: 1:50-1:100
Molecular Weight	-

Background

Cooperates with LY96 and CD14 to mediate the innate immune response to bacterial lipopolysaccharide (LPS).


Product Name: Toll-Like Receptor 4 Rabbit Polyclonal Antibody
Catalog #: APRab00443



Research Area

Immunology

Image Data



Immunohistochemistry analysis of paraffin-embedded mouse Spleen Tissue using Toll-Like Receptor 4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.