

Product Name: CD232 Rabbit Polyclonal Antibody
Catalog #: APRab00568



Summary

Production Name	CD232 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	PLXNC1
Alternative Names	PLXNC1; VESPR; Plexin-C1; Virus-encoded semaphorin protein receptor; CD232
Gene ID	10154
SwissProt ID	O60486.

Application

Dilution Ratio	WB: 1:500-1:1000 ELISA: 1:10000
Molecular Weight	Calculated MW: 176 kDa; Observed MW: 176 kDa

Background

Receptor for SEMA7A, for smallpox semaphorin A39R, vaccinia virus semaphorin A39R and for herpesvirus Sema protein.


Product Name: CD232 Rabbit Polyclonal Antibody
Catalog #: APRab00568




Research Area

Cancer

Image Data



Western blot analysis of CD232 in 3T3 lysates using CD232 antibody.



Western blot analysis of CD232 in 3T3 lysates using CD232 antibody.

Note

For research use only.