

Product Name: Epigen Rabbit Polyclonal Antibody
Catalog #: APRab00572



Summary

Production Name	Epigen Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	EPGN
Alternative Names	EPGN; Epigen; Epithelial mitogen; EPG
Gene ID	255324
SwissProt ID	Q6UW88.

Application

Dilution Ratio	WB: 1:500-1:1000 ELISA: 1:10000
Molecular Weight	Calculated MW: 17 kDa; Observed MW: 17 kDa

Background

Promotes the growth of epithelial cells.


Product Name: Epigen Rabbit Polyclonal Antibody
Catalog #: APRab00572




Research Area

Signal Transduction

Image Data



Western blot analysis of Epigen in 293T, 823 lysates using Epigen antibody.



Western blot analysis of Epigen in 293T lysates using Epigen antibody.

Note

For research use only.