

Product Name: Integrin alpha 3 Rabbit Polyclonal Antibody
Catalog #: APRab00576

Summary

Production Name	Integrin alpha 3 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,IHC-P,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	ITGA3
Alternative Names	ITGA3; MSK18; Integrin alpha-3; CD49 antigen-like family member C; FRP-2; Galactoprotein B3; GAPB3; VLA-3 subunit alpha; CD49c
Gene ID	3675
SwissProt ID	P26006.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 ELISA: 1:10000
Molecular Weight	Calculated MW: 117 kDa; Observed MW: 117 kDa

Product Name: Integrin alpha 3 Rabbit Polyclonal Antibody
Catalog #: APRab00576




Background

Integrin alpha-3/beta-1 is a receptor for fibronectin, laminin, collagen, epiligrin, thrombospondin and CSPG4.


Research Area

Signal Transduction

Image Data



Western blot analysis of Integrin alpha 3 in hela lysates using Integrin alpha 3 antibody.



Immunohistochemistry analysis of paraffin-embedded mouse spleen using Integrin alpha 3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.