

Product Name: HGS Rabbit Polyclonal Antibody
Catalog #: APRab00612



Summary

Production Name	HGS Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	HGS
Alternative Names	HGS; HRS; Hepatocyte growth factor-regulated tyrosine kinase substrate; Hrs; Protein pp110
Gene ID	9146
SwissProt ID	O14964.

Application

Dilution Ratio	WB: 1:500-1:1000 ELISA: 1:10000
Molecular Weight	Calculated MW: 86 kDa; Observed MW: 86 kDa

Background

Involved in intracellular signal transduction mediated by cytokines and growth factors.


Product Name: HGS Rabbit Polyclonal Antibody
Catalog #: APRab00612



Research Area

Cell Biology

Image Data



Western blot analysis of HGS in mouse eye, mouse kidney, mouse testis lysates using HGS antibody.

Note

For research use only.