

Product Name: TCEAL6 Rabbit Polyclonal Antibody
Catalog #: APRab00636



Summary

Production Name	TCEAL6 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	TCEAL6
Alternative Names	TCEAL6; Transcription elongation factor A protein-like 6; TCEA-like protein 6; Transcription elongation factor S-II protein-like 6
Gene ID	158931.0
SwissProt ID	Q6IPX3.

Application

Dilution Ratio	WB: 1:500-1:1000 ELISA: 1:10000
Molecular Weight	Calculated MW: 22 kDa; Observed MW: 32 kDa

Background

May be involved in transcriptional regulation.


Product Name: TCEAL6 Rabbit Polyclonal Antibody
Catalog #: APRab00636



Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of TCEAL6 in CoLo lysates using TCEAL6 antibody.

Note

For research use only.