

Product Name: PI3 Kinase C2 gamma Rabbit Polyclonal Antibody
Catalog #: APRab00639



Summary

Production Name	PI3 Kinase C2 gamma Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	PIK3C2G PIK3C2G; Phosphatidylinositol 4-phosphate 3-kinase C2 domain-containing subunit
Alternative Names	gamma; PI3K-C2-gamma; PtdIns-3-kinase C2 subunit gamma; Phosphoinositide 3-kinase-C2-gamma
Gene ID	5288
SwissProt ID	O75747.

Application

Dilution Ratio	WB: 1:500-1:1000 ELISA: 1:10000
Molecular Weight	Calculated MW: 166 kDa; Observed MW: 160 kDa

Product Name: PI3 Kinase C2 gamma Rabbit Polyclonal Antibody
Catalog #: APRab00639




Background

Generates phosphatidylinositol 3-phosphate (PtdIns3P) and phosphatidylinositol 3,4-bisphosphate (PtdIns(3,4)P₂) that act as second messengers (By similarity).

Research Area

Cell Biology

Image Data



Western blot analysis of PI3 Kinase C2 gamma in A549 lysates using PI3 Kinase C2 gamma antibody.

Note

For research use only.