

Product Name: Collagen XXVIII alpha 1 Rabbit Polyclonal Antibody
Catalog #: APRab00645

Summary

Production Name	Collagen XXVIII alpha 1 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	COL28A1
Alternative Names	Collagen alpha-1(XXVIII) chain
Gene ID	340267
SwissProt ID	Q2UY09.

Application

Dilution Ratio	WB: 1:500-1:1000 ELISA: 1:10000
Molecular Weight	Calculated MW: 117 kDa; Observed MW: 117 kDa

Background


May act as a cell-binding protein.

**Product Name: Collagen XXVIII alpha 1 Rabbit
Polyclonal Antibody
Catalog #: APRab00645**

Research Area

Cell Biology

Image Data



Western blot analysis of Collagen XXVIII alpha 1 in 3T3 HEPG2 lysates using Collagen XXVIII alpha 1 antibody.

Note

For research use only.