

Product Name: MonoMethyl-Histone H3 (Lys36) Rabbit Polyclonal Antibody
Catalog #: APRab00688

Summary

Production Name	MonoMethyl-Histone H3 (Lys36) Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,ICC/IF
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Methylated
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification	Affinity Chromatography

Immunogen

Gene Name	H3C1
Alternative Names	H3K36me; H3/j; H3C1; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FJ; H3C10; H3C11; HIST1H3J
Gene ID	8350
SwissProt ID	P68431.

Application

Dilution Ratio	WB: 1:500-1:1000 IF: 1:50-1:200
Molecular Weight	Calculated MW: 15 kDa; Observed MW: 15 kDa

**Product Name: MonoMethyl-Histone H3 (Lys36) Rabbit
Polyclonal Antibody**
Catalog #: APRab00688


Background

H3 Core component of nucleosome. Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessibility to the cellular machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of Histone H3 (mono methyl K36) in (1) HeLa lysates; (2) NIH/3T3 lysates using MonoMethyl-Histone H3 (Lys36) antibody.

Note

For research use only.