

Product Name: B3GALT2 Rabbit Polyclonal Antibody
Catalog #: APRab00815



Summary

Production Name	B3GALT2 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	IHC-F,IHC-P,ICC/IF,ELISA
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	B3GALT2
Alternative Names	B3GALT2; Beta-1; 3-galactosyltransferase 2; Beta-1; 3-GalTase 2; Beta3Gal-T2; Beta3GalT2; UDP-galactose:2-acetamido-2-deoxy-D-glucose 3beta-galactosyltransferase 2
Gene ID	8707
SwissProt ID	O43825.

Application

Dilution Ratio	IHC: 1:50-1:100 IF: 1:50-1:200 ELISA: 1:10000
Molecular Weight	-

Background

Product Name: B3GALT2 Rabbit Polyclonal Antibody
Catalog #: APRab00815




Beta-1,3-galactosyltransferase that transfers galactose from UDP-galactose to substrates with a terminal beta-N-acetylglucosamine (beta-GlcNAc) residue.

Research Area

Cell Biology

Image Data



Immunohistochemistry analysis of paraffin-embedded Human heart using B3GALT2 antibody. High-pressure and temperature Tris-EDTA pH 8.0 was used for antigen retrieval. Sample with blocking peptide on the right.

Note

For research use only.