

Product Name: TriMethyl-Histone H3 (Lys14) Rabbit Polyclonal Antibody
Catalog #: APRab00850



Summary

Production Name	TriMethyl-Histone H3 (Lys14) Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Methylated
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	H3-4
Alternative Names	H3K14me3; H3 histone; HIST1H3A; Histone cluster 1; H3a
Gene ID	8290
SwissProt ID	Q16695.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 16 kDa; Observed MW: 16 kDa

Background

Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in

**Product Name: TriMethyl-Histone H3 (Lys14) Rabbit
Polyclonal Antibody**
Catalog #: AP Rab00850




eukaryotes. Nucleosomes consist of approximately 146 bp of DNA wrapped around a histone octamer composed of pairs of each of the four core histones (H2A, H2B, H3, and H4). The chromatin fiber is further compacted through the interaction of a linker histone, H1, with the DNA between the nucleosomes to form higher order chromatin structures. This gene is intronless and encodes a member of the histone H3 family. Transcripts from this gene lack polyA tails; instead, they contain a palindromic termination element. This gene is located separately from the other H3 genes that are in the histone gene cluster on chromosome 6p22-p21.3.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of TriMethyl-Histone H3 (Lys14) in various cell lines lysates using TriMethyl-Histone H3K14 antibody.

Note

For research use only.