

Product Name: Acetyl-Lysine Rabbit Polyclonal Antibody
Catalog #: APRab03406

Summary

Production Name	Acetyl-Lysine Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Species independent

Performance

Conjugation	Unconjugated
Modification	Acetylated
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	-
Alternative Names	-
Gene ID	-
SwissProt ID	-.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Refer to figures

Background


Acetylation of lysine residues within proteins has emerged as an important mechanism used by cells to overcome this repression.

Product Name: Acetyl-Lysine Rabbit Polyclonal Antibody
Catalog #: APRab03406

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of Acetyl-Lysine in HeLa lysates using Pan Acetyl-Lysine antibody.

Note

For research use only.