

Product Name: Rrn3 (phospho Ser649) Rabbit Polyclonal Antibody
Catalog #: APRab05390

Summary

Production Name	Rrn3 (phospho Ser649) Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,IHC-P,IF-P,IF-F,ICC/IF,ELISA
Reactivity	Human,Rat,Mouse

Performance

Conjugation	Unconjugated
Modification	Phospho Antibody
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	RRN3
Alternative Names	RRN3; TIFIA; RNA polymerase I-specific transcription initiation factor RRN3; Transcription initiation factor IA; TIF-IA
Gene ID	54700.0
SwissProt ID	Q9NYV6.The antiserum was produced against synthesized peptide derived from human TIF-IA around the phosphorylation site of Ser649. AA range:602-651

Application

Dilution Ratio	WB 1:500-1:2000, IHC-P 1:100-1:300, ELISA 1:10000, IF-P/IF-F/ICC/IF 1:50-200
Molecular Weight	70kDa


Product Name: Rrn3 (phospho Ser649) Rabbit Polyclonal Antibody
Catalog #: APRab05390

Background


function:Required for efficient transcription initiation by RNA polymerase I.,similarity:Belongs to the RRN3 family.,function:Required for efficient transcription initiation by RNA polymerase I.,similarity:Belongs to the RRN3 family.,

Research Area

Image Data




Enzyme-Linked Immunosorbent Assay (Phospho-ELISA) for Immunogen Phosphopeptide (Phospho-left) and Non-Phosphopeptide (Phospho-right) , using TIF-IA (Phospho-Ser649) Antibody



Immunohistochemistry analysis of paraffin-embedded human brain, using TIF-IA (Phospho-Ser649) Antibody. The picture on the right is blocked with the phospho peptide.

Product Name: Rrn3 (phospho Ser649) Rabbit Polyclonal Antibody
Catalog #: APRab05390



Western blot analysis of lysates from Jurkat cells treated with starved 24h, using TIF-IA (Phospho-Ser649) Antibody. The lane on the right is blocked with the phospho peptide.

Note

For research use only.