

## Summary

<b>Production Name</b>	Calregulin Rabbit Polyclonal Antibody
<b>Description</b>	Rabbit Polyclonal Antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,FC,IHC-P,IF-P,IF-F,ICC/IF,ELISA
<b>Reactivity</b>	Human,Mouse,Rat,Monkey

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
<b>Purification</b>	Affinity purification

## Immunogen

<b>Gene Name</b>	CALR
<b>Alternative Names</b>	CALR; CRTC; Calreticulin; CRP55; Calregulin; Endoplasmic reticulum resident protein 60; ERp60; HACBP; grp60
<b>Gene ID</b>	811.0
<b>SwissProt ID</b>	P27797.The antiserum was produced against synthesized peptide derived from human CALR. AA range:21-70

## Application

<b>Dilution Ratio</b>	WB 1:500-2000, Flow Cyt 1:50-200, IHC-P 1:100-500, IF-P/IF-F/ICC/IF/ICC 1:100-500, ELISA 1:5000-20000
<b>Molecular Weight</b>	48kDa


## Background

Calreticulin is a multifunctional protein that acts as a major  $\text{Ca}^{2+}$ -binding (storage) protein in the lumen of the endoplasmic reticulum. It is also found in the nucleus, suggesting that it may have a role in transcription regulation. Calreticulin binds to the synthetic peptide KLGFFKR, which is almost identical to an amino acid sequence in the DNA-binding domain of the superfamily of nuclear receptors. Calreticulin binds to antibodies in certain sera of systemic lupus and Sjogren patients which contain anti-Ro/SSA antibodies, it is highly conserved among species, and it is located in the endoplasmic and sarcoplasmic reticulum where it may bind calcium. The amino terminus of calreticulin interacts with the DNA-binding domain of the glucocorticoid receptor and prevents the receptor from binding to its specific glucocorticoid response element. Calreticulin can inhibit the binding of androgen receptor to its caution: Was originally (PubMed:2332496) thought to be the 52 kDa Ro autoantigen., domain: Associates with PDIA3 through the tip of the extended arm formed by the P-domain., domain: Can be divided into a N-terminal globular domain, a proline-rich P-domain forming an elongated arm-like structure and a C-terminal acidic domain. The P-domain binds one molecule of calcium with high affinity, whereas the acidic C-domain binds multiple calcium ions with low affinity., domain: The interaction with glycans occurs through a binding site in the globular lectin domain., domain: The zinc binding sites are localized to the N-domain., function: Molecular calcium binding chaperone promoting folding, oligomeric assembly and quality control in the ER via the calreticulin/calnexin cycle. This lectin interacts transiently with almost all of the monoglucosylated glycoproteins that are synthesized in the ER. Interacts with the DNA-binding domain of NR3C1 and mediates its nuclear export., mass spectrometry: PubMed:11149926, online information: Calreticulin, online information: Calreticulin entry, similarity: Belongs to the calreticulin family., subcellular location: Also found in cell surface (T cells), cytosol and extracellular matrix. Associated with the lytic granules in the cytolytic T-lymphocytes., subunit: Monomer. Component of an EIF2 complex at least composed of CUGBP1, CALR, CALR3, EIF2S1, EIF2S2, HSP90B1 and HSPA5. Interacts with PDIA3/ERp57 and with NR3C1.,

## Research Area

Antigen processing and presentation;


## Image Data




**Product Name: Calregulin Rabbit Polyclonal Antibody**  
**Catalog #: APRab07874**




Immunofluorescence analysis of NIH/3T3 cells, using CALR Antibody. The picture on the right is blocked with the synthesized peptide.



Immunohistochemistry analysis of paraffin-embedded human heart tissue, using CALR Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from COS7 cells, using CALR Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using Calregulin Polyclonal Antibody

**Note**

For research use only.