

Product Name: GAAP Rabbit Polyclonal Antibody
Catalog #: APRab11222



Summary

Production Name	GAAP Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,IHC-P,IF-P,IF-F,ICC/IF,ELISA
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	TMBIM4
Alternative Names	TMBIM4; GAAP; LFG4; CGI-119; Protein lifeguard 4; Golgi anti-apoptotic protein; Protein S1R; Transmembrane BAX inhibitor motif-containing protein 4; Z-protein
Gene ID	51643.0
SwissProt ID	Q9HC24. The antiserum was produced against synthesized peptide derived from human TMBIM4. AA range:121-170

Application

Dilution Ratio	WB 1:500-1:2000, IHC-P 1:100-1:300, ELISA 1:5000, IF-P/IF-F/ICC/IF 1:50-200
Molecular Weight	23kDa

Background


Product Name: GAAP Rabbit Polyclonal Antibody
Catalog #: APRab11222



similarity:Belongs to the B11 family.,similarity:Belongs to the B11 family.,

Research Area

Image Data



Western blot analysis of lysates from HepG2 cells, using TMBIM4 Antibody. The lane on the right is blocked with the synthesized peptide.

Note

For research use only.