

Product Name: GOT1B Rabbit Polyclonal Antibody
Catalog #: APRab11584



Summary

Production Name	GOT1B Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	GOLT1B GOT1A CGI-141 HDCMA39P UNQ432/PRO793
Alternative Names	
Gene ID	51026.0
SwissProt ID	Q9Y3E0. Synthesized peptide derived from human protein . at AA range: 70-150

Application

Dilution Ratio	WB 1:500-2000, ELISA 1:5000-20000
Molecular Weight	15kDa

Background

function:May be involved in fusion of ER-derived transport vesicles with the Golgi complex.,similarity:Belongs to the GOT1 family.,function:May be involved in fusion of ER-derived transport vesicles with the Golgi complex.,similarity:Belongs to the


Product Name: GOT1B Rabbit Polyclonal Antibody
Catalog #: APRab11584



GOT1 family.,

Research Area

Image Data



Western blot analysis of lysates from U2OS cells, primary antibody was diluted at 1:1000, 4°over night

Note

For research use only.