

Product Name: GRP1 Rabbit Polyclonal Antibody
Catalog #: APRab11791



Summary

Production Name	GRP1 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse,Rat,Monkey

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	CYTH3 CYTH3; ARNO3; GRP1; PSCD3; Cytohesin-3; ARF nucleotide-binding site opener 3;
Alternative Names	Protein ARNO3; General receptor of phosphoinositides 1; Grp1; PH; SEC7 and coiled-coil domain-containing protein 3
Gene ID	9265.0
SwissProt ID	O43739.The antiserum was produced against synthesized peptide derived from human GRP1. AA range:351-400

Application


Dilution Ratio	WB 1:500-1:2000, ELISA 1:10000.Not yet tested in other applications.
Molecular Weight	46kDa

Background


This gene encodes a member of the PSCD (pleckstrin homology, Sec7 and coiled-coil domains) family. PSCD family members have identical structural organization that consists of an N-terminal coiled-coil motif, a central Sec7 domain, and a C-terminal pleckstrin homology (PH) domain. The coiled-coil motif is involved in homodimerization, the Sec7 domain contains guanine-nucleotide exchange protein (GEP) activity, and the PH domain interacts with phospholipids and is responsible for association of PSCDs with membranes. Members of this family appear to mediate the regulation of protein sorting and membrane trafficking. This encoded protein is involved in the control of Golgi structure and function, and it may have a physiological role in regulating ADP-ribosylation factor protein 6 (ARF) functions, in addition to acting on ARF1. [provided by RefSeq, Jul 2008],function:Promotes guanine-nucleotide exchange on ARF1. Promotes the activation of ARF through replacement of GDP with GTP.,similarity:Contains 1 PH domain.,similarity:Contains 1 SEC7 domain.,subunit:Binds via its PH domain to the inositol head group of phosphatidylinositol 3,4,5-triphosphate with high affinity. Interacts with GRASP.,tissue specificity:Almost absent from liver, thymus and peripheral blood lymphocytes.,

Research Area

Image Data



Western blot analysis of lysate from COS7 cells, using GRP1 antibody.



Product Name: GRP1 Rabbit Polyclonal Antibody
Catalog #: APRab11791



Western Blot analysis of various cells using GRP1 Polyclonal Antibody

Note

For research use only.