

**Product Name: Histone deacetylase 6 Rabbit Polyclonal Antibody**  
**Catalog #: APRab12048**

---

## Summary

<b>Production Name</b>	Histone deacetylase 6 Rabbit Polyclonal Antibody
<b>Description</b>	Rabbit Polyclonal Antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IHC-P,IF-P,IF-F,ICC/IF,ELISA
<b>Reactivity</b>	Human,Mouse

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
<b>Purification</b>	Affinity purification

## Immunogen

<b>Gene Name</b>	HDAC6
<b>Alternative Names</b>	HDAC6; KIAA0901; JM21; Histone deacetylase 6; HD6
<b>Gene ID</b>	10013.0
<b>SwissProt ID</b>	Q9UBN7.The antiserum was produced against synthesized peptide derived from human HDAC6. AA range:1166-1215

## Application

<b>Dilution Ratio</b>	WB 1:500-1:2000, IHC-P 1:100-1:300, IF-P/IF-F/ICC/IF 1:200-1:1000, ELISA 1:20000.Not yet tested in other applications.
<b>Molecular Weight</b>	135kDa

**Product Name: Histone deacetylase 6 Rabbit Polyclonal Antibody**  
**Catalog #: AP Rab12048**

---


## Background

Histones play a critical role in transcriptional regulation, cell cycle progression, and developmental events. Histone acetylation/deacetylation alters chromosome structure and affects transcription factor access to DNA. The protein encoded by this gene belongs to class II of the histone deacetylase/acuc/apha family. It contains an internal duplication of two catalytic domains which appear to function independently of each other. This protein possesses histone deacetylase activity and represses transcription. [provided by RefSeq, Jul 2008], catalytic activity: Hydrolysis of an N(6)-acetyl-lysine residue of a histone to yield a deacetylated histone., function: Responsible for the deacetylation of lysine residues on the N-terminal part of the core histones (H2A, H2B, H3 and H4). Histone deacetylation gives a tag for epigenetic repression and plays an important role in transcriptional regulation, cell cycle progression and developmental events. Histone deacetylases act via the formation of large multiprotein complexes (By similarity). Plays a central role in microtubule-dependent cell motility via deacetylation of tubulin., PTM: Sumoylated in vitro., PTM: Ubiquitinated. Its polyubiquitination however does not lead to its degradation., similarity: Belongs to the histone deacetylase family. Type 2 subfamily., similarity: Contains 1 UBP-type zinc finger., subcellular location: It is mainly cytoplasmic, where it is associated with microtubules., subunit: Interacts with CBFA2T3, HDAC11 and SIRT2. Interacts with F-actin. Interacts with BBIP10.,

## Research Area


Protein\_Acetylation

## Image Data




Immunohistochemistry analysis of paraffin-embedded human colon carcinoma tissue, using HDAC6 Antibody. The picture on the right is blocked with the synthesized peptide.


**Product Name: Histone deacetylase 6 Rabbit Polyclonal Antibody**  
**Catalog #: APRab12048**



Western blot analysis of lysates from HeLa cells, using HDAC6 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using Histone deacetylase 6 Polyclonal Antibody diluted at 1: 2000



Western Blot analysis of MCF7 cells using Histone deacetylase 6 Polyclonal Antibody diluted at 1: 2000

## Note

For research use only.