

Product Name: ISK7 Rabbit Polyclonal Antibody
Catalog #: APRab12771



Summary

Production Name	ISK7 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	SPINK7 ECG2 UNQ745/PRO1474
Alternative Names	Serine protease inhibitor Kazal-type 7 (Esophagus cancer-related gene 2 protein) (ECRG-2)
Gene ID	84651.0
SwissProt ID	P58062.

Application

Dilution Ratio	WB 1:500-2000
Molecular Weight	


Background

Product Name: ISK7 Rabbit Polyclonal Antibody
Catalog #: APRab12771



Research Area

Image Data



Western Blot analysis of hela 293T using primary antibody at 1:1000 dilution 4°C, overnight. Secondary antibody was diluted at 1:10000 25°C, 1.5hours

Note

For research use only.