

Summary

Production Name	M-RIP Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,ELISA,IHC-P
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	MPRIP
Alternative Names	MPRIP; KIAA0864; MRIP; RHOIP3; Myosin phosphatase Rho-interacting protein; M-RIP; Rho-interacting protein 3; RIP3; p116Rip
Gene ID	23164/23164
SwissProt ID	Q6WCQ1.Synthesized peptide derived from the Internal region of human M-RIP.

Application

Dilution Ratio	WB 1:500-2000, IHC-P 1:50-300, ELISA 2000-20000
Molecular Weight	115kDa

Background


Product Name: M-RIP Rabbit Polyclonal Antibody
Catalog #: APRab14095




function: Targets myosin phosphatase to the actin cytoskeleton. Required for the regulation of the actin cytoskeleton by RhoA and ROCK1. Depletion leads to an increased number of stress fibers in smooth muscle cells through stabilization of actin fibers by phosphorylated myosin. Overexpression of MRIP as well as its F-actin-binding region leads to disassembly of stress fibers in neuronal cells., similarity: Contains 1 PH domain., similarity: Contains 2 PH domains., subcellular location: Colocalizes with F-actin., subunit: Binds F-actin through its N-terminus (By similarity). Binds RHOA, PPP1R12A/MBS and PPP1R12C/MBS85 through adjacent coiled coil domains., function: Targets myosin phosphatase to the actin cytoskeleton. Required for the regulation of the actin cytoskeleton by RhoA and ROCK1. Depletion leads to an increased number of stress fibers in smooth muscle cells through stabilization of actin fibers by phosphorylated myosin. Overexpression of MRIP as well as its F-actin-binding region leads to disassembly of stress fibers in neuronal cells., similarity: Contains 1 PH domain., similarity: Contains 2 PH domains., subcellular location: Colocalizes with F-actin., subunit: Binds F-actin through its N-terminus (By similarity). Binds RHOA, PPP1R12A/MBS and PPP1R12C/MBS85 through adjacent coiled coil domains.,

Research Area

Image Data



Western Blot analysis of various cells using M-RIP Polyclonal Antibody



Western Blot analysis of 293 cells using M-RIP Polyclonal Antibody

Product Name: M-RIP Rabbit Polyclonal Antibody
Catalog #: APRab14095



Note

For research use only.