

Product Name: NY-REN-50 Rabbit Polyclonal Antibody
Catalog #: APRab15005



Summary

Production Name	NY-REN-50 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	IHC-P,IF-P,IF-F,ICC/IF,ELISA
Reactivity	Human,Rat,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	ZNF839
Alternative Names	ZNF839; C14orf131; Zinc finger protein 839; Renal carcinoma antigen NY-REN-50
Gene ID	55778.0
SwissProt ID	A8K0R7. Synthesized peptide derived from the Internal region of human NY-REN-50.

Application

Dilution Ratio	IHC-P 1:100-1:300, ELISA 1:40000, IF-P/IF-F/ICC/IF 1:50-200
Molecular Weight	

Background


similarity:Contains 1 C2H2-type zinc finger.,similarity:Contains 1 C2H2-type zinc finger.,

Product Name: NY-REN-50 Rabbit Polyclonal Antibody
Catalog #: APRab15005



Research Area

Image Data



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200 (4° overnight) . 2, Tris-EDTA, pH 9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200 (room temperature, 30min) .

Note

For research use only.