

Product Name: Olfactory receptor 11G2 Rabbit Polyclonal Antibody
Catalog #: APRab15154

Summary

Production Name	Olfactory receptor 11G2 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,IF-P,IF-F,ICC/IF,ELISA
Reactivity	Human,Rat,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	OR11G2
Alternative Names	OR11G2; Olfactory receptor 11G2; Olfactory receptor OR14-34
Gene ID	390439.0
SwissProt ID	Q8NGC1.The antiserum was produced against synthesized peptide derived from human OR11G2. AA range:241-290

Application

Dilution Ratio	WB 1:500-1:2000, IF-P/IF-F/ICC/IF 1:200-1:1000, ELISA 1:40000.Not yet tested in other applications.
Molecular Weight	35kDa

Product Name: Olfactory receptor 11G2 Rabbit Polyclonal Antibody
Catalog #: APRab15154


Background

Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008],caution:It is uncertain whether Met-1 or Met-35 is the initiator.,function:Odorant receptor ,similarity:Belongs to the G-protein coupled receptor 1 family.,


Research Area

Olfactory transduction;

Image Data




Immunofluorescence analysis of A549 cells, using OR11G2 Antibody. The picture on the right is blocked with the synthesized




Western blot analysis of lysates from HUVEC, Jurkat, and MCF-7 cells, using OR11G2 Antibody. The lane on the right is

**Product Name: Olfactory receptor 11G2 Rabbit
Polyclonal Antibody
Catalog #: APRab15154**

blocked with the synthesized peptide.



Western blot analysis of the lysates from HeLa cells using OR11G2 antibody.



Western Blot analysis of various cells using Olfactory receptor 11G2 Polyclonal Antibody diluted at 1: 1000

Note

For research use only.