

Product Name: Olfactory receptor 56B1 Rabbit Polyclonal Antibody
Catalog #: APRab15272

Summary

Production Name	Olfactory receptor 56B1 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	IF-P,IF-F,ICC/IF,ELISA
Reactivity	Human,Rat,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	OR56B1
Alternative Names	OR56B1; OR56B1P; Olfactory receptor 56B1; Olfactory receptor OR11-65
Gene ID	387748.0
SwissProt ID	Q8NGI3. The antiserum was produced against synthesized peptide derived from human OR56B1. AA range:161-210

Application

Dilution Ratio	IF-P/IF-F/ICC/IF 1:200-1:1000, ELISA 1:10000. Not yet tested in other applications.
Molecular Weight	

Background


**Product Name: Olfactory receptor 56B1 Rabbit
Polyclonal Antibody
Catalog #: APRab15272**

Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008],caution:It is uncertain whether Met-1 or Met-4 is the initiator.,function:Odorant receptor ,similarity:Belongs to the G-protein coupled receptor 1 family,.

Research Area

Olfactory transduction;

Image Data



Immunofluorescence analysis of HUVEC cells, using OR56B1 Antibody. The picture on the right is blocked with the synthesized peptide.

Note

For research use only.