

## Summary

<b>Production Name</b>	PKC $\gamma$ Rabbit Polyclonal Antibody
<b>Description</b>	Rabbit Polyclonal Antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IHC-P,IF-P,IF-F,ICC/IF,ELISA
<b>Reactivity</b>	Human,Mouse,Rat

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
<b>Purification</b>	Affinity purification

## Immunogen

<b>Gene Name</b>	PRKCG
<b>Alternative Names</b>	PRKCG; PKCG; Protein kinase C gamma type; PKC-gamma
<b>Gene ID</b>	5582.0
<b>SwissProt ID</b>	P05129.The antiserum was produced against synthesized peptide derived from the Internal region of human PRKCG. AA range:521-570

## Application

<b>Dilution Ratio</b>	WB 1:500-1:2000, IHC-P 1:100-300, ELISA 1:20000, IF-P/IF-F/ICC/IF 1:50-200
<b>Molecular Weight</b>	78kDa

## Background

**Product Name: PKC  $\gamma$  Rabbit Polyclonal Antibody**  
**Catalog #: APRab16201**

---




Protein kinase C (PKC) is a family of serine- and threonine-specific protein kinases that can be activated by calcium and second messenger diacylglycerol. PKC family members phosphorylate a wide variety of protein targets and are known to be involved in diverse cellular signaling pathways. PKC also serve as major receptors for phorbol esters, a class of tumor promoters. Each member of the PKC family has a specific expression profile and is believed to play distinct roles in cells. The protein encoded by this gene is one of the PKC family members. This protein kinase is expressed solely in the brain and spinal cord and its localization is restricted to neurons. It has been demonstrated that several neuronal functions, including long term potentiation (LTP) and long term depression (LTD), specifically require this kinase. Knockout studies in mice also suggest that this kinase may be involved in neuromodulatory activity: ATP + a protein = ADP + a phosphoprotein., cofactor: Binds 3 calcium ions per subunit. The ions are bound to the C2 domain., disease: Defects in PRKCG are the cause of spinocerebellar ataxia type 14 (SCA14) [MIM:605361]. Spinocerebellar ataxia is a clinically and genetically heterogeneous group of cerebellar disorders. Patients show progressive incoordination of gait and often poor coordination of hands, speech and eye movements, due to degeneration of the cerebellum with variable involvement of the brainstem and spinal cord. SCA14 is an autosomal dominant cerebellar ataxia (ADCA), function: PKC is activated by diacylglycerol which in turn phosphorylates a range of cellular proteins. PKC also serves as the receptor for phorbol esters, a class of tumor promoters., function: This is a calcium-activated, phospholipid-dependent, serine- and threonine-specific enzyme., online information: Retina International's Scientific Newsletter, similarity: Belongs to the protein kinase superfamily. AGC Ser/Thr protein kinase family. PKC subfamily., similarity: Contains 1 AGC-kinase C-terminal domain., similarity: Contains 1 C2 domain., similarity: Contains 1 protein kinase domain., similarity: Contains 2 phorbol-ester/DAG-type zinc fingers., subunit: Interacts with CDCP1.,

## Research Area


MAPK\_ERK\_Growth; MAPK\_G\_Protein; ErbB\_HER; Calcium; Phosphatidylinositol signaling system; Vascular smooth muscle contraction; WNT; WNT-T CELLVEGF; Focal adhesion; Tight junction; Gap junction; Natural killer cell mediated cytotoxicity; Fc gamma R-mediated phagocytosis; Leukocyte transendothelial migration; Long-term potentiation; Long-term depression; Melanogenesis; Aldosterone-regulated sodium reabsorption; Vibrio cholerae infection; Pathways in cancer; Glioma; Non-small cell lung cancer;

## Image Data


**Product Name: PKC  $\gamma$  Rabbit Polyclonal Antibody**  
**Catalog #: APRab16201**



Western blot analysis of lysate from A549 cells, using PRKCG Antibody.



Western Blot analysis of A549, THP-1, K562, mouse brain, mouse skeletal muscle, mouse testis cells using PKC  $\gamma$  Polyclonal Antibody.. Secondary antibody was diluted at 1:20000



Immunohistochemical analysis of paraffin-embedded rat-brain, antibody was diluted at 1:100

## Note

For research use only.