

Product Name: Ribosomal Protein L7L1 Rabbit Polyclonal Antibody
Catalog #: APRab17170

Summary

Production Name	Ribosomal Protein L7L1 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,IHC-P,IF-P,IF-F,ICC/IF,ELISA
Reactivity	Human,Rat,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	RPL7L1
Alternative Names	RPL7L1; 60S ribosomal protein L7-like 1
Gene ID	285855.0
SwissProt ID	Q6DKI1. The antiserum was produced against synthesized peptide derived from human RPL7L1. AA range:41-90

Application

Dilution Ratio	WB 1:500-1:2000, IHC-P 1:100-1:300, ELISA 1:40000, IF-P/IF-F/ICC/IF 1:50-200
Molecular Weight	29kDa


Background

**Product Name: Ribosomal Protein L7L1 Rabbit
Polyclonal Antibody
Catalog #: APRab17170**

similarity:Belongs to the ribosomal protein L30P family.,similarity:Belongs to the ribosomal protein L30P family.,

Research Area

Image Data



Note

For research use only.