

Summary

Production Name	SR-5A Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,IF-P,IF-F,ICC/IF,ELISA
Reactivity	Human,Rat,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	HTR5A
Alternative Names	HTR5A; 5-hydroxytryptamine receptor 5A; 5-HT-5; 5-HT-5A; 5-HT5A; Serotonin receptor 5A
Gene ID	3361.0
SwissProt ID	P47898.The antiserum was produced against synthesized peptide derived from human HTR5A. AA range:15-64

Application

Dilution Ratio	WB 1:500-1:2000, IF-P/IF-F/ICC/IF 1:200-1:1000, ELISA 1:10000.Not yet tested in other applications.
Molecular Weight	36kDa


Background

The neurotransmitter serotonin (5-hydroxytryptamine, 5-HT) has been implicated in a wide range of psychiatric conditions and also has vasoconstrictive and vasodilatory effects. The gene described in this record is a member of 5-hydroxytryptamine (serotonin) receptor family and encodes a multi-pass membrane protein that functions as a receptor for 5-hydroxytryptamine and couples to G-proteins. This protein has been shown to function in part through the regulation of intracellular Ca²⁺ mobilization. [provided by RefSeq, Jul 2008],function:This is one of the several different receptors for 5-hydroxytryptamine (serotonin), a biogenic hormone that functions as a neurotransmitter, a hormone, and a mitogen. The activity of this receptor is mediated by G proteins.,similarity:Belongs to the G-protein coupled receptor 1 family,.


Research Area

Calcium;Neuroactive ligand-receptor interaction;

Image Data




Immunofluorescence analysis of LOVO cells, using HTR5A Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from Jurkat cells, using HTR5A Antibody. The lane on the right is blocked with the synthesized peptide.

Product Name: SR-5A Rabbit Polyclonal Antibody
Catalog #: APRab18253



Western blot analysis of the lysates from HepG2 cells using HTR5A antibody.

Note

For research use only.