

Product Name: T3JAM Rabbit Polyclonal Antibody
Catalog #: APRab18581



Summary

Production Name	T3JAM Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,IHC-P,IF-P,IF-F,ICC/IF,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	TRAF3IP3
Alternative Names	TRAF3IP3; T3JAM; TRAF3-interacting JNK-activating modulator; TRAF3-interacting protein 3
Gene ID	80342.0
SwissProt ID	Q9Y228.The antiserum was produced against synthesized peptide derived from human T3JAM. AA range:251-300

Application

Dilution Ratio	WB 1:500-1:2000, IHC-P 1:100-1:300, IF-P/IF-F/ICC/IF 1:200-1:1000, ELISA 1:20000.Not yet tested in other applications.
Molecular Weight	64kDa

Product Name: T3JAM Rabbit Polyclonal Antibody
Catalog #: APRab18581




Background


The gene encodes a protein that mediates cell growth by modulating the c-Jun N-terminal kinase signal transduction pathway. The encoded protein may also interact with a large multi-protein assembly containing the phosphatase 2A catalytic subunit. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2013],function:May function as an adapter molecule that regulates TRAF3-mediated JNK activation.,subunit: Binds to the isoleucine zipper of TRAF3 via its coiled-coil domain.,

Research Area

Image Data




Immunofluorescence analysis of HeLa cells, using T3JAM Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from Jurkat cells, using T3JAM Antibody. The lane on the right is blocked with the synthesized peptide.

Product Name: T3JAM Rabbit Polyclonal Antibody
Catalog #: APRab18581



Western Blot analysis of various cells using T3JAM Polyclonal Antibody

Note

For research use only.